 The EU directive 93/42/EEC was applied in the design and production of this medical device.

**NAKANISHI INC.** 

[www.nsk-inc.com](http://www.nsk-inc.com)

700 Shimohinata Kanuma-shi  
Tochigi 322-8666,  
Japan

**NSK Europe GmbH** 

[www.nsk-europe.de](http://www.nsk-europe.de)

Eily-Beinhorn-Strasse 8  
65760 Eschborn,  
Germany

**NSK United Kingdom Ltd**

[www.nsk-uk.com](http://www.nsk-uk.com)

Office 5, Gateway1000,  
Arlington Business Park, Whittle Way,  
Stevenage, SG1 2FP, UK

**NSK America Corp**

[www.nskdental.com](http://www.nskdental.com)

1800 Global Parkway  
Hoffman Estates, IL 60192,  
USA

**NSK France SAS**

[www.nsk.fr](http://www.nsk.fr)

32 rue de Lisbonne  
75008 Paris,  
France

**NSK Dental Spain SA**

[www.nsk-spain.es](http://www.nsk-spain.es)

C/ M6dena, 43 El Soho-Eur6polis  
28232 Las Rozas, Madrid,  
Spain

**NSK Oceania Pty Ltd**

[www.nskoceania.com.au](http://www.nskoceania.com.au)

Unit 22, 198-222 Young St.  
Waterloo, Sydney,  
NSW 2017, Australia

**NSK Middle East**

[www.nsk-inc.com](http://www.nsk-inc.com)

Room 6EA-701, 7th Floor, East Wing No.6  
Dubai Airport Free Zone,  
PO Box 54316 Dubai, UAE

**NSK Asia Pte Ltd**

[www.nsk-inc.com](http://www.nsk-inc.com)

1 Maritime Square,  
#09-33 HarbourFront Centre,  
Singapore 099253

Specifications are subject to change without notice.

'11.08.02 ©

OPERATION MANUAL

**NSK**

 0197

Please read this Operation Manual carefully  
before use, and file for future reference

**Ti-Max**

Electric Motor

**M40-LED  
M40B-LED**

**MADE IN JAPAN**

OM-E0463E 002

Powerful Partners®

## Intended to use

This motor is designed only for dental clinic use by qualified personnel.

### CAUTION

- When operating the motor always consider the safety of the patient.
- Do not allow any impact and do not drop the motor.
- Check the vibration, noise and overheating before use and should the motor function abnormally, stop using immediately and contact dealer for repair.
- Refer to the motor cord manual for connecting it to the motor.
- Do not connect / disconnect the attachment until the drive motor has completely stopped.
- Do not exceed the speed recommended by the handpiece manufacture.
- Do not lubricate the motor.
- Do not wipe with, or clean or immerse in, high acid water or sterilizing solutions.

**CAUTION**

- Do not sterilize a high temperature including the Autoclave sterilization.
- Perform periodical maintenance checks.
- If the motor has not been used for a long period, check for noise, rotation, and overheating before use.

## 1. Connecting / Disconnecting the Motor Cord

### (1) Connecting

- 1) Align the motor and the pipe hole of the motor cord.
- 2) Insert the Motor.
- 3) Fasten the Motor Cord Nut securely.

### (2) Disconnecting

Unscrew and detach the Motor Cord Nut, and pull out the Motor Cord Connector.


Fig. 1

**NOTICE for M40B-LED**

Remove Motor Cord Adapter before connecting.

## 2. Connecting / Disconnecting the Handpiece

### (1) Connecting

- 1) Insert the E-type handpiece in the Motor Insert.
- 2) Turn the handpiece until it clicks. (Positioning mark is aligned.)


Fig. 2

- (2) Disconnecting  
Simply pull out the handpiece from the motor.

**CAUTION**

Care should be taken when using any handpiece attachment which requires lubrication. Following lubrication stand the handpiece attachment vertically until the excess solution drains. Wipe the handpiece then attach it to the motor. If oil enters to the motor, it may cause malfunction.

### 3. Cleaning

After each use, clean the surface of the motor with brush etc. (Do not use the metal brush), and wipe it with the cotton moistened with disinfecting alcohol.

### 4. Replace the O-rings at the Motor Insert

- (1) Remove the worn O-rings with a pointed tool like a needle.  
(2) Insert new O-rings the O-ring grooves.

\* Optional O-ring (Blue) : Order No. 0313 084070  
Optional O-ring : Order No. 0312 074080

**CAUTION**

There are 4 O-Rings on the Motor Insert. Blue one is thinner than other 3 pieces. Make sure to order it.


Fig. 3


Fig. 4

## Specification

Model	M40-LED	M40B-LED
Speed	60-40,000min <sup>-1</sup> (rpm)	
Dimensions	Ø21 x L103mm	Ø21 x L91.75mm
Hose Connection	Bien-Air® 4VLM and 4VR 400	
Optics	LED	
Voltage	AC/DC 3.3V ± 0.05V	
Activate Voltage	2.5 - 4.0 V	
Consumption Current	0.35 Tpy. (3.3V)	

For Lighting power, follow the EN60601-1, Voltage have to be match SELV (Safety Extra Low Voltage)  
The guide light and illuminating radiation area may be influenced by conditions of light volume and color.

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