



Dentalone

The dentist's office Mobile cart



NAKANISHI INC.

www.nsk-inc.com

700 Shimohinata
Kanuma-shi
Tochigi 322-8666,
Japan

TEL : +81 (0)289-64-3380
FAX : +81 (0)289-62-5636

EC REP

NSK Europe GmbH

www.nsk-europe.de

Elly-Beinhorn-Strasse 8
65760 Eschborn,
Germany

TEL : +49 (0)6196 77606-0
FAX : +49 (0)6196 77606-29

NSK United Kingdom Ltd

www.nsk-uk.com

Office 5, Gateway1000,
Arlington Business Park,
Whittle Way,
Stevenage, SG1 2FP, UK

TEL : +44 (0)1438 310670
FAX : +44 (0)1438 310671

NSK Dental LLC

www.nskdental.com

1800 Global Parkway
Hoffman Estates, IL 60192,
USA

TEL : +1 888 675 1675
FAX : +1 800 838 9328

NSK Oceania Pty Ltd

www.nskoceania.com.au

Unit 22, 198-222 Young St.
Waterloo, Sydney,
NSW 2017, Australia

TEL : +61 (0)2 8306 3000
FAX : +61 (0)2 8306 3030

NSK Middle East

www.nsk-inc.com

Room 6EA-701, 7th Floor,
East Wing No.6
Dubai Airport Free Zone,
PO Box 54316 Dubai, UAE

TEL : +971 (0)4-204-5507
FAX : +971 (0)4-204-5508

NSK Asia Pte Ltd

www.nsk-inc.com

1 Maritime Square,
#09-33 HarbourFront Centre,
Singapore 099253

TEL : +65-6376-5018
FAX : +65-6376-0233

NSK France

www.nsk.fr

32, rue de Lisbonne
75008 Paris
France

TEL : +33(0)1 567 959 80
FAX : +33(0)1 567 959 81

NSK Dental Spain SA

www.nsk-spain.es

C/ Módena, 43 El Soho - Európolis
28232 Las Rozas, Madrid - España
TEL : +34-91-626-61-28
FAX : +34-91-626-61-32



Army & Hospitals



Humanitarian



Home Care



Schools

EXTENDED FUNCTIONALITIES

ADVANTAGES

- Compact
- Simple to use
- Mobile
- Plug and Play system
- Easy to carry
- User-friendly
- Extensible to Cart function
- Large range functionalities

Dentalone is the dentist's office mobile cart easily expandable from «carry» position to a comfortable, user-friendly dental treatment cart ready to use anywhere, anytime.

Dentalone offers a large range of instruments and extended functionalities, thus practitioners can realise dental treatments outside any dentist's office:

- Instruments carrier with automatic selection of instruments
- Illuminated LCD control panel
- 1 Brushless NLX nano micromotor with LED Light
- 1 Integrated Multi Function Ultrasonic Scaler Varios 170 with LED Light
- 1 Syringe with air/water/spray with LED Light
- 1 Cannula for saliva

DENTIST'S OFFICE MOBILE CART

Plug and play



User-friendly



FEATURES

CONTENT

- 1 Carrying Case
- 1 Water tank for spray
- 1 Tank for vacuum
- Built in compressor
- Built in suction pump
- 1 Variable Foot Control
- 1 Extendable Tablet

SPECIFICATIONS

Dimensions (HxWxD)

- Case with trolley (unpulled trolley handle): 545 x 351 x 405 mm
- Case with trolley (pulled trolley handle): 1075 x 351 x 405 mm
- Without case: 496 x 350 x 405 mm
- Working height: 766 mm
- Weight with case: 18,5 kg
- Power Supply: 100 - 230 V~ 50/60 Hz

Compact



Expandable to Cart Function

496mm



766mm



ACCESSORIES

MODEL	ORDER CODE	PRODUCTS
Saliva complete with hose and connector	57064	
3F syringe with light	57066	
Screwing cap suction, assembled	57068	
Suction bottle complete	57072	
Screwing cap water, assembled	57069	
Water tank complete	57073	

MODEL	ORDER CODE	PRODUCTS
Cover assembled	57082	
Socket assembled	57088	
NLX motor cord	37074	
Cord for Scaler Varios 170	37075	
Foot control Z1082001X	37100	
Sterilization case	Z1035001	
Tip Wrench	Z221076	



Extended Functionalities



Easy to transport