



No	Type	No	Type	No	Type
<b>gmi</b>		<b>Implantvision</b>		<b>Medimecca</b>	
avantgard RP 3.75-4.25	55066 58	Internal Hex 3.7-5.7	55034 26	PT SLA Chaorum	55047 39
avantgard WP 4.75-5.75	55075 67				
frontier RP,WP	55051 43	<b>Intoss Anchor</b>		<b>Megagen</b>	
phoenix NP 3.3	55038 30	Internal Hex 3.7-5.7	55088 80	Anyridge	55051 43
phoenix Std 3.3 - Wide 5.0	55011 1			Anyone Internal	55037 29
Monolith 3.0	55043 35	<b>Intra-lock</b>		Mini	55064 56
		Intrahex 3.4 - 4.5	55051 43		
<b>GP Implants</b>		Gold & Blue 3.4 - 6.0	55051 43	<b>Microdent</b>	
SSI Spiral Shape implant	55051 43	InDex 3.4 - 3.75	55015 5	Universal	55011 1
		InDex 4.0 - 4.3	55020 10	Microdent 2.8 - 3.25	55068 60
<b>gt medical</b>		FlatTop 3.75	55012 2	Microdent 3.3 - 5.0	55011 1
HXE	55011 1			Microdent 5.5 - 8.0	55013 3
HXE RP	55012 2	<b>IRES</b>		Genius 3.0	55025 16
HXE WP	55013 3	External Hex Narrow 3.3	55012 2	Genius 3.5 - 5.0	55051 43
HXI 4.1	55037 29	External Hex Regular 3.7 - 5.2	55011 1	Ektos 3.2	55023 14
002 HXI 4.7	55049 41	Internal Hex Narrow 3.3	55087 79	Ektos 3.7 - 8.0	55051 43
OCT 4.3	55014 4	Internal Hex Regular 3.7 - 5.2	55049 41	Trylogic 3.5 - 4.2	55051 43
		Internal Octagon Regular 3.7 - 5.2	55014 4	Trylogic 5.0	55075 67
<b>Hahn</b>				<b>MIS</b>	
Tapered 3.0	55058 50	<b>J Dental Care</b>		Conical connection Narrow (C1, V3)	55081 73
Tapered 3.5	55046 38	JDEvolution 3.7-6.0	55067 59	Conical connection Standard (C1, V3)	55087 79
Tapered 4.3	55046 38	JD Evolution S 3.25	55044 36	Conical connection Wide (C1)	55041 33
Tapered 5.0	55037 29	JD ICON ultra, S 2.75	55064 56	Seven and M4 Narrow	55027 19
Tapered 7.0	55048 40	Conical Abutment	55080 72	Seven and M4 Standard	55049 41
				Seven Wide 5.0	55051 43
<b>Hiossen</b>		<b>Keystone</b>			
HS System	55028 20	Prima / Genesis Narrow	55049 41		
ET System	55037 29	Prima / Genesis Regular	55067 59	<b>Neo Biotech</b>	
ET 3.2	55047 39	Prima / Genesis Wide	55041 33	IS II 4.0	55020 10
ET mini	55015 5	Restore Narrow	55079 71	IS III 4.0	55020 10
		Restore Regular	55011 1		
<b>Hi-Tec</b>		Restore Wide	55013 3	<b>Neodent</b>	
Tapered Self Thread NTI 3.3, SLT 3,75	55034 26			Drive CM	55025 16
Tapered Self Thread GTI 4.2, WTI 5.0	55051 43	<b>Klockner</b>		Titamax WS Cortical	55056 48
Spark TPI	55034 26	VEGA Micro 3.0	55074 66	GM	55066 58
Logic Plus	55034 26	VEGA Narrow 3.5	55065 57		
Logic LGI 3.0-3.5	55046 38	VEGA Regular 4.0, 4.5	55034 26	<b>Neoss</b>	
Logic LGI 4.3-5.0	55037 29	KL Narrow 3.3	55012 2	SP (∅ 3.5 implants and wider)	55004 18
		KL Regular 3.7, 4.2	55011 1	NP (∅ 3.25)	55027 19
<b>IDI Evolution</b>		KL Wide 4.7	55013 3	Access Abutment	55060 52
Standard implant 3.75-4.0	55004 18	SK2	55038 30		
		Essential Cone	55067 59	<b>Nobel Biocare</b>	
		S3M 3.1	55079 71	Brånemark System Ext hex RP	55011 1
<b>IDI Implants</b>				Brånemark System Ext hex NP	55012 2
ID Cam ST 4.2	55010 13	<b>Konmet</b>		Brånemark System Ext hex WP	55013 3
ID All	55010 13	Classic 4.0, 4.8	55014 4	Tri-lobe Connection NP	55022 12
		Cylindrical 3.5-4.0	55065 57	Tri-lobe Connection RP, WP	55010 13
<b>ids</b>		Cylindrical 4.5-5.5	55057 49	Tri-lobe Connection 6.0	55032 24
Reflect Aspire	55081 73	Short Cylindrical 5.0-6.0	55078 70	Conical Connection 3.0	55058 50
Reflect Certus	55019 9	Conical 4.4/3.4	55065 57	Conical Connection Narrow	55046 38
Reflect Tapered	55034 26	Conical 4.8/3.8, 5.5/4.5	55057 49	Conical Connection Regular	55037 29
Reflect Recover	55020 10			Conical Connection Wide 5.5	55048 40
Reflect Rapid	55020 10	<b>Lasak</b>		MUA NP & RP	55080 72
		BioniQ QR	55066 58	On1 NP, RP, WP	55042 34
<b>Implance</b>		BioniQ QN	55021 11		
Tissue level	55014 4			<b>Noris Medical</b>	
Bone level Regular	55037 29	<b>Leader</b>		Tuff 3.3 - 6.0	55034 26
Bone level 3.3	55065 57	Internal Hex, 3.3 - 5.5 mm	55034 26		
		Cono-In	55067 59	<b>Novem</b>	
<b>Implant Direct</b>				3.4	55039 31
Legacy 1/2/3/4 3.0/3.2	55012 2			3.8	55016 6
Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051 43	<b>Maxillent</b>		4.2, 5.0	55057 49
RePlant 3.5	55022 12	iRaise	55049 41		
RePlant 4.3 - 5.0	55010 13			<b>NucleOss</b>	
RePlant 6.0	55032 24	<b>Medentika</b>		Tpure Bone Level	55004 18
ReActive 3.5 Platform	55022 12	Microcone NI	55021 11	T6 Bonel Level	55074 66
ReActive 4.3 - 5.0 Platform	55010 13	Microcone RI	55081 73		
Interactive 3.0 Platform	55046 38	Quattrocone RI	55081 73	<b>OCO Biomedical</b>	
Interactive 3.4 Platform	55037 29	Quattrocone 30 AI	55081 73	Engage 3.25	55027 19
Screwindirect 5.0 PF 3.0 - 5.7	55070 62	Procone 3.8	55031 23	Engage 4.0-6.0	55051 43
		Procone 4.3	55032 24	ERI 3.25	55069 61
		Procone 5.0		ERI 4.0-5.0	55014 4
<b>Implant One</b>				TSI 3.25	55068 60
100 Series Std thread	55029 21	<b>Medentis</b>		TSI 4.0-5.0	55011 1
200 Series Std thread	55029 21	ICX Premium 3.75 - 4.8	55062 54		
300 Series Std thread (blue)	55029 21	ICX Plus 3.45	55039 31	<b>Osstem</b>	
400 Series Std thread (purple)	55035 27	ICX Narrow	55058 50	US Regular	55011 1
500 series Std thread (green)	55035 27			US Mini	55012 2
300 Series Wide thread (blue)	55029 21	<b>Medical Instinct</b>		US Wide	55013 3
400 Series Wide thread (purple)	55035 27	Hex 3.0	55071 63	US Wide PS	55011 1
500 Series Wide thread (green)	55035 27	Hex 3.4, 4.0, 5.0	55052 44	SS Regular, Wide	55028 20
		Cone 3.0, 3.4, 4.0	55071 63	TS Regular, Reg Short, Diam 6 & 7	55037 29
<b>Implant Swiss</b>		Cone 5.0	55072 64	TS Mini 3.5	55015 5
Tissue level	55014 4			TS Mini 3.0	55047 39
Bone level Regular	55037 29				
Bone level 3.3	55065 57				

No	Type	No	Type	No	Type
<b>Osteo Plus</b>		<b>Simply Implants</b>		<b>Thommen</b>	
Iridium PF 3.8 & 4.2	55079 71	Internal Hex	55051 43	SPI 3.5	55021 11
She 3.45, 3.75	55011 1			SPI 4.0	55065 57
Shi 3.75, 4.2	55049 41	<b>SIN</b>		SPI 4.5-6.0	55066 58
Shi 5.0	55051 43	Strong SW CM	55023 14		
		Strong SW Int Hex	55035 27	<b>ticare</b>	
<b>OT-Medical</b>		Unitite CM 2.9	55084 76	Osseus Mini, Standard	55011 1
OT-F2 3.4, 3.8	55053 45	Unitite CM 3.5	55077 69	Osseus Maxi	55013 3
OT-F2 4.1, 5.0	55073 65	Unitite CM Compact 4.0	55085 77	Osseus Maxi ps	55011 1
		Tryon 3.25 - 5.0	55011 1	Inhex Mini	55021 11
				Inhex Standard	55046 38
<b>Overmed</b>				Inhex Maxi	55017 7
Easyline Easydip	55020 10	<b>Southern Implants</b>		Osseus Quattro Standard	55011 1
		Tri-Nex 4.3-5.0	55010 13	Inhex Quattro Standard	55046 38
<b>Oxy Implant</b>		Tri-Nex 3.5	55022 12		
Internal	55020 10	External Hex	55011 1		
Kone Conical narrow	55021 11	MAX Ext Hex	55011 1	<b>TRI Dental implants</b>	
Kone Conical Regular	55051 43	Co-Axis Ext Hex	55011 1	Tissue level Octa	55014 4
<b>PSK Line Regular</b>	<b>55051 43</b>	Piccolo Ext Hex	55050 42	Tri-Vent, all	55034 26
<b>External Mini</b>	<b>55079 71</b>	IT-Connection	55014 4	Narrow	55025 16
<b>External Regular</b>	<b>55011 1</b>	Deep Conical 3.0	55064 56		
<b>External Maxi</b>	<b>55011 1</b>	Deep Conical 3.5-4.0	55046 38	<b>Trinon</b>	
		Deep Conical Tapered 5.0	55017 7	Q2-Implant	55011 1
				QK-Implant	55014 4
<b>Paltop</b>				QZA 3.5 Platform	55049 41
Standard 3.75-5.0	55034 26	<b>Spiral Tech</b>		QZA 4.5, 5.7 Platform	55051 43
Narrow	55027 19	ESi, Hex connection	55051 43		
PCA Conical	55066 58	Ultimate, Hex connection	55051 43		
		ESi Conical 3.5	55044 36	<b>Vulkan</b>	
		ESi Conical 4.3, 5.0	55037 29	Internal Hex	55034 26
<b>Phibo</b>				External Hex NP	55012 2
Aurea EVO NP	55027 19			External Hex RP	55011 1
Aurea EVO RP, WP	55051 43	<b>Straumann</b>		Conical con MiP	55058 50
TSA S3	55074 66	Tissue level, Regular & Wide Neck	55014 4	Conical con NP	55046 38
TSA S4, S5	55034 26	Tissue level, Narrow Neck Crossfit (NNC)	55076 68		
TSH S3, S4, S5	55011 1	Bone level / Bone Level Tapered Narrow	55061 53		
		Bone level / Bone Level Tapered Regular	55062 54	<b>Winsix</b>	
		Bone level / Bone Level Tapered SC 2.9	55064 56	Internal Hex	55020 10
<b>Prodent Italia</b>		Bone Level X (BLX) all platform	55027 19		
Twinner & Prime, Conomet, Collar: Orange	55027 19	Screw retained Abutment NP, RP	55080 72	<b>Xive</b>	
Twinner & Prime, Conomet, Collar: Fuchsia	55051 43			S-Plus 3.0	55030 22
<b>Twinner &amp; Prime, Conomet, Collar: Green</b>	<b>55067 59</b>			S-Plus 3.4	55061 53
Twinner & Prime, Conomet, Collar: Yellow	55020 10	<b>Surgikor</b>		S-Plus 3.8-5.5	55062 54
Prime, Conomet, Collar: Blue	55020 10	Fixation 3.0	55058 50	TG Plus 3.4-4.5	55018 8
		Fixation 3.5	55066 58		
		Fixation 4.3	55037 29		
<b>Quest</b>		Solution 3.5, 4.0	55066 58	<b>Zeramex</b>	
Hylock NP	55030 22	Solution 4.5, 6.0	55037 29	XT Implant 4.2, 5.5	55086 78
Hylock RP	55066 58	Immediate, all	55049 41		
		Versatile, all	55049 41	<b>Zimmer</b>	
<b>Resista</b>				Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049 41
WIS 3.0-5.0	55049 41			Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051 43
		<b>Sweden &amp; Martina</b>		Tapered Swiss Platf 3.8	55070 62
<b>Rex Implants</b>		Syra, Syra SP 3.9-5.0	55011 1	Tapered Swiss Platf 4.8	55014 4
External hex TL 1.8	55060 52	Outlink <sup>2</sup> 3.75, 4.10, 5.0	55011 1	Spline Twist Platf 4 & 5	55038 30
Internal hex TL 2.9	55034 26	Shelta 3.8-5.0	55025 16	Eztetic NP 3.1	55065 57
		Premium Straight & SP 3.3-3.8	55025 16		
		Kohno Straight & SP 3.3-3.8	55025 16		
<b>Ritter</b>		Premium Straight & SP 4.25-6.0	55004 18	<b>Zuga</b>	
Spiral SB/LA Narrow PF 2.9	55081 73	Kohno Straight & SP 4.25-6.0	55004 18	All	55027 19
Spiral SB/LA Standard PF 3.75	55049 41	Prama 3.8 - 5.0	55079 71		
		<b>DAT Connection</b>	<b>55035 27</b>		
<b>RS Implant</b>					
3.5 - 6.0	55020 10				
		<b>Tag Dental</b>			
<b>Schutz Dental</b>		Axis	55051 43		
Impla cylindrical	55052 44	massif	55051 43		
IMPLA micro Retention Cone 4.2	55052 44				
		<b>Tatum</b>			
<b>SGS Implant</b>		Integrity 3.7 - 5.0	55015 5		
P1, P7 Platform	55049 41	Integrity 6.0 - 8.0	55040 32		
P1D, P7D	55049 41	T implant	55063 55		
P7N	55059 51				
		<b>TBR</b>			
<b>SIC</b>		Infinity PF 3.5, PF 4.0	55004 18		
Ace 3.4 - 5.0	55052 44	M Implant	55044 36		
Vantage 3.0	55026 17				
Vantage 4.2-5.2	55054 46				